

resiDNA_workflows


Residual DNA qPCR Analysis — Colab Notebook and Local Streamlit App


Project Organization

```
resiDNA_workflows/
|
├─ app/
|   └─ app.py                # Phase II Streamlit app – local web UI, mirrors
|   ↳ the notebook
|   └─ audit_log.py          # opt-in audit log (RESIDNA_AUDIT_LOG) –
|   ↳ app.py-only, see Data Privacy below
|   └─ requirements.txt      # dependencies for the Streamlit app
├─ .streamlit/
|   └─ config.toml           # light theme + wide layout (must stay at repo
|   ↳ root)
|
├─ Notebooks/
|   └─ 01_Phase I.ipynb      # Phase I notebook – import through PDF report
|   ↳ generation
|
├─ Data/
|   └─ experiment_001.xlsx    # primary dataset the notebook was built against
|   └─ sample_data_001.xls   # real-world dataset – validated
|   └─ sample_data_002.xls   # real-world dataset – validated
|
├─ Scripts/
|   └─ qpcr.py               # shared classification/suitability logic –
|   ↳ imported by both the notebook and app.py
|   └─ plotting.py           # shared style_table() – imported by both the
|   ↳ notebook and app.py
├─ Figures/                  # reserved – future shared modules
├─ Results/                  # generated PDF reports
├─ Logs/                     # audit.log, only if RESIDNA_AUDIT_LOG is enabled
|   ↳ – gitignored, not committed
├─ README.md
├─ IT_Deployment_Guide.md    # IT handoff: on-prem / Azure + SharePoint
|   ↳ hosting paths (Option A decided)
├─ IT_Setup_Guide_OptionA.md # from-scratch runbook IT executes: commands,
|   ↳ service files, config templates
```

Application Architecture (In Progress)

Design of a local-first architecture, transitioning to IT-hosted per IT_Deployment_Guide.md.





Development Road Map

Phase I - Colab Scientific Prototype (.ipynb) >Prove that the qPCR analysis logic works correctly.

1. [x] Load QuantStudio file 2. [x] QuantStudio Results Parser 3. [x] Data Classification Engine 4. [x] System Suitability Engine 5. [x] Sample Result Processing 6. [x] Create Prototype Graphs 7. [x] Validate against Data/sample_data_001.xls and Data/sample_data_002.xls 8. [x] PDF Report generation

Phase II - Convert Prototype into a Local Application (<http://localhost:8501>) >Create a simple website you run on your own computer using Local Streamlit App. 1. [x] Upload Page 2. [x] Backend Processing 3. [x] Results Dashboard

Phase III - Add Criteria Management >Stop changing code every time acceptance criteria change. >Update thresholds directly in the notebook's #@param cell and mirror them into app.py — no separate Settings UI or database needed.

Phase IV - Improve Reporting and Visualization — Not needed >Current interactive Plotly charts and the PDF report already meet the bar for a professional analytical tool; dropped from the roadmap.

Phase V - Hand Off to Company IT — Next >Move from prototype to maintainable, IT-hosted software. >Pin dependencies, organize into a professional VS Code project, and transition hosting to company-controlled infrastructure per IT_Deployment_Guide.md (on-prem/local server or Azure + SharePoint embed) with Entra ID SSO.

Phase VI - (Maybe) Upgrade to Full Web Application >Only worth revisiting if IT's hosted deployment outgrows what Streamlit can offer (e.g. multi-user concurrency, custom auth flows beyond OIDC). >Replace Streamlit with React Frontend and FastAPI Backend

Acceptance Criteria

STD (Standard Curve) | Check | Formula | Pass | |—|—|—| | R^2 | instrument-computed | ≥ 0.99 | | Ct %CV (triplicate) | $(Ct_SD / Ct_Mean) \times 100 \leq 20\%$ | | Back-calculation bias | $(Back-Calc Qty - Nominal Qty) / Nominal Qty \times 100 \leq 25\%$ | | Amplification efficiency | $(10^{(-1/Slope)} - 1) \times 100 \leq 90\%-110\%$ |

ERC & PC (Controls) | Check | Formula | Pass | |—|—|—| | Quantity %CV (triplicate) | $(Qty_SD / Qty_Mean) \times 100 \leq 20\%$ | | PC %Recovery | $(Dilution-Adjusted Qty / Reference-Adjusted Qty) \times 100 \leq 80\%-125\%$ | | ERC %Recovery | same formula | 50%–150% |

Sample LOQ / STD Range — replaces the old ICH-sigma LOQ formula. From the standard curve's own back-calculated values, take STD1 (highest nominal quantity) and STD6 (lowest): | Value | Source | Used for | |—|—|—| | LOQ (Quantity) | STD6 Back-Calc Mean | Sample STD Range / “below LOQ” floor | | LOQ (Ct) | STD6 Ct Mean | NTC / NEC pass check | | STD Range | [STD6 Back-Calc Mean, STD1 Back-Calc Mean] | Sample STD Range check |

NTC / NEC — both pass per well if Ct is Undetermined or $Ct \geq LQ_Ct$ (STD6's own measured Ct Mean above).

Sample Suitability — evaluated in this order, since each step feeds the next:

1. Triplicate single-outlier exclusion: if a dilution's full 3-well Quantity %CV fails 25%, `resolve_sample_replicates()` (`Scripts/qpcr.py`) drops the outlier well and re-checks the remaining 2; if they pass, that pair's mean becomes the dilution's Quantity %CV / Quantity Mean / Total DNA / DNA per Protein (“Replicates Used: 2”) instead of an outright fail. Only triggers on an already-failing triplicate — passing ones keep the instrument's reported figures as-is. Every check below uses this resolved Quantity Mean.
2. Quantity %CV: $(Qty_SD / Qty_Mean) \times 100 \leq 25\%$.
3. STD Range: the resolved Quantity Mean (pre-dilution-adjustment) must fall within the STD Range above; below it is Out of Range and is “below LOQ”, above it is also Out of Range.
4. Dilutional Linearity: $\%Bias = (x_1 - x_2) / x_2 \times 100 \leq 20\%$, where x_2 is the Dilution Adjusted quantity of the sample's own dilution with the lowest combined Ct %CV + Quantity %CV among dilutions that already pass both step 2 and step 3, and x_1 is each dilution's Dilution Adjusted quantity. A sample with no dilution passing steps 2–3 has no valid reference, so every one of its dilutions gets N/A.
5. Combined Suitability: Pass only if steps 2–4 all pass; otherwise Fail (binary — no N/A, since an uninterpretable dilution isn't a soft middle ground). Only Pass dilutions get averaged into the final result.

Two tables: Final Sample Results by Dilution and Final Sample Results — Averaged per Sample (every sample is still reported, even one with zero passing dilutions — averaged across its passing dilutions if any exist, otherwise across all of them as a flagged fallback). In the averaged table: red highlight = the sample failed acceptance criteria (takes priority), green highlight = passed and DNA per Protein is below SAMPLE_DNA_PER_PROTEIN_LIMIT (15 ng/mg), gray + – LOQ suffix on Total DNA / DNA per Protein = the averaged Quantity Mean is below the LOQ (low confidence, footnoted). Sample # (e.g. S1) groups dilutions into the averaged table but isn't itself displayed in the by-dilution table — Sample (e.g. S1 D1) already carries that information.

Next Step — `compute_sample_status()` (`Scripts/qpcr.py`) turns the by-dilution Suitability breakdown into a next-action label per sample: | Next Step | Meaning | |—|—| | Reportable | At least one dilution cleanly passes all of steps 2–4 above (same condition as Sample Passed). || LOQ – Spike Test | Quantity %CV and Dilutional Linearity pass (Linearity N/A counts as non-blocking), but STD Range fails because the dilution's Quantity Mean is below STD6 — below the assay's LOQ. || ULQ – Adjust Dilution | Same as above, but STD Range fails because Quantity Mean is above STD1 — above the curve's highest calibrated point. || Retest | Everything else: STD Range passes but Quantity %CV or Linearity fails, or multiple checks fail together. |

Each dilution is classified individually first, then a sample with multiple dilutions reports whichever status ranks best across them (Reportable > LOQ – Spike Test / ULQ – Adjust Dilution > Retest) — the same “any dilution passing is enough” rule already used for Sample Passed.

Criteria Verification

- Output Verification: all thresholds live in one Colab-form-editable cell (#@param) in the notebook, right after Data Classification. Values marked “instrument-computed” above are read forward from the QuantStudio export as-is (not recalculated) — the notebook's Output Verification cell checks they

survive parsing unchanged. The one exception is Dilutional Linearity's % Bias, which is genuinely computed rather than read from the instrument.

- Manual Golden-Reference Check: a one-time sanity check for that computed % Bias value. Hand-calculate (or Excel-calculate) the expected % Bias for a few known dilutions against the reference-dilution formula, enter them into the notebook's GOLDEN_LINEARITY_BIAS dict, and re-run the cell — it flags any dilution where the notebook's computed value doesn't match your expected value within a ± 0.01 percentage-point tolerance (GOLDEN_BIAS_TOLERANCE).

PDF Report

Both the notebook and the Streamlit app generate the same formatted PDF summarizing a run — Run Info, System Suitability, Sample Suitability, Final Sample Results, and a Signatures block. In the notebook, this is the final ## PDF Report section, using the #@param fields in the Report Info cell; it writes to Results/Sample_Analysis_Report_<run_no>.pdf and requires reportlab (auto-installs if missing). In the app, it's built the same way but in-memory and offered as a browser download instead.

- The logo is currently a placeholder box in both — swap in a real logo file.

Streamlit App (DEMO)

Demo deployment: <https://residna.streamlit.app> — hosted on Streamlit Community Cloud, no installation needed, demo use only. Access is restricted to authorized company viewers (see Data Privacy & Security before uploading real company data). See below to run it locally instead.

app.py is the Phase II local web app: upload a QuantStudio file and get the same System Suitability, Sample Suitability, and Final Sample Results sections as the notebook, without touching code, plus: - Sample ID / Sample Name inputs — one row per Base Sample actually detected in the upload (simpler than the notebook's 8 fixed #@param slots, since the app can generate inputs from real data). - Standard Curve graphs (QC-controls overlay and per-sample overlay) and Amplification Curves (ΔR_n vs Cycle) — all interactive Plotly charts. - A Download PDF Report button producing the same report described above.

Only the sections that map to the PDF report are shown — raw file preview, parse stats, and the classification table run silently in the background rather than cluttering the page.

Data Privacy & Security (DEMO Streamlit App)

In transit - The Community Cloud link is served over HTTPS by default — traffic between the browser and Streamlit's servers is encrypted.

In the app - app.py processes uploaded files entirely in memory (st.file_uploader → pandas parsing → Plotly charts / PDF generation). By default, nothing is written to disk or a database on the server — no Results/ writes, no SQLite. Once a session ends or the app restarts, the uploaded file and its analysis are gone. - The PDF report is streamed straight to the browser as a download; it isn't retained server-side. - Exception — audit log: app/audit_log.py appends one record per generated PDF report (timestamp, uploaded filename, a SHA-256 hash of it, assay/run info, submitter/reviewer name, pass/fail — never the analytical results themselves) to Logs/audit_log, but only if the RESIDNA_AUDIT_LOG environment variable is set. It's unset (off) by default and on this demo deployment, so the “nothing persisted” behavior above still holds here — it's only enabled on the company-hosted deployment, where compliance requires it (see IT_Deployment_Guide.md).

Who can reach the app - Access to the live link is restricted via Streamlit Community Cloud's viewer authentication: in the app's dashboard under Settings → Sharing, it is set to an allowlist of company email addresses — only people who sign in with one of those emails can open it. - Residual risk worth knowing: even restricted, uploaded data still passes through Streamlit's shared Community Cloud infrastructure, which is intended for demos/public apps rather than as a substitute for company-controlled hosting. For the most sensitive runs, run the app locally instead (see below) so data never leaves your machine or network.

For moving to a fully company-hosted deployment (on-prem/local server or Azure + SharePoint embed), see IT_Deployment_Guide.md.

Running the Notebook Locally

The notebook auto-detects Colab vs. local execution (IN_COLAB flag) — outside Colab it skips `git clone/pip install klib/google.colab` imports and resolves paths relative to the repo instead of `/content/`.

```
python -m venv .venv
.venv\Scripts\pip install pandas numpy openpyxl xlrd klib ipykernel plotly
↪ reportlab
.venv\Scripts\python -m ipykernel install --user --name resiDNA-venv
↪ --display-name "resiDNA (.venv)"
```

Then select the resiDNA (.venv) kernel in VS Code / Jupyter and run all cells.

Running the Streamlit App Locally

Run these from the repo root — `.streamlit/config.toml` only applies when Streamlit is launched from the directory that contains it.

```
python -m venv .venv
.venv\Scripts\pip install -r app\requirements.txt
.venv\Scripts\streamlit run app\app.py
```

Opens at `http://localhost:8501` — only reachable while that command is still running in a terminal; closing it (or restarting your computer) stops the app until you run it again.

The audit log (see Data Privacy above) is off by default locally too. To try it: set `RESIDNA_AUDIT_LOG=1` before launching (`$env:RESIDNA_AUDIT_LOG="1"` in PowerShell), generate a report, then check `Logs/audit.log`.

Running Tests

```
python -m venv .venv
.venv\Scripts\pip install -r app\requirements.txt -r tests\requirements.txt
.venv\Scripts\pytest tests -v
```

Three suites, all under `tests/`: - `test_qpcr.py` — unit tests for every pure function in `Scripts/qpcr.py`, independent of any data file. - `test_app_golden.py` — runs the real `app/app.py` (headlessly, via a fake streamlit module in `tests/streamlit_stub.py`) against each file in `Data/` and checks `final_results/reportable_results` against known-good values, so a future change that silently shifts a Sample Suitability result gets caught automatically. `Data/*.xls/.xlsx` are real (scrubbed) company data and aren't committed to the repo (see Data Privacy above and git history), so these tests skip, not fail, wherever those files aren't present — a fresh clone or CI run will show them skipped, which is expected. - `test_audit_log.py` — unit tests for `app/audit_log.py`, including that it's a no-op unless `RESIDNA_AUDIT_LOG` is set and that raw analytical data never ends up in a log record.

Backlog — Recommended Cleanup

- ☐ Notebook outputs (styled tables, the `klib` plot) are committed inline, which makes diffs large. Consider `nbstripout` or committing a rendered HTML export alongside a stripped notebook.
- ☐ `Scripts/qpcr.py` / `Scripts/plotting.py` now hold the shared, self-contained logic (`style_table`, `classify_sample`, `recovery_bounds`, `combined_suitability`, the LOD sigma-Ct formula, the single-outlier triplicate resolution (`resolve_sample_replicates`), the Dilutional Linearity %Bias

calculation, the Sample LOQ/STD Range derivation, STD Range suitability, the averaged-results fallback aggregation, and the per-sample Next Step classification) that both the notebook and `app.py` import. Deliberately left duplicated: the parsing/cleaning steps and suitability-table assembly, since the notebook keeps those inline, cell-by-cell, on purpose — that step-by-step visibility is Phase I's whole point (see Road Map). Worth revisiting only if that tradeoff stops being worth it.